


Health Monitoring System



Motivation

- Develop an efficient system that can monitor the patients' health parameters like temperature and pulse rate .
- System that can effectively deliver the data to a patient monitoring system where it is stored permanently in a database.
- Investigating the potential of Arduino and Raspberrypi in collecting, sending multiple sensed parameters in real-time.



Sensing

- This stage is responsible for reading patient's temperature and pulse rate.
- The sensors used are
 - Temperature sensor- MAX90614 IR Sensor
 - Pulse oximeter sensor- MAX30100 Pulse Oximeter
- Microcontrollers used are
 - NodeMCU
 - RaspberryPi

Networking

- Data sensed by the sensors is captured by the microcontrollers.
- This data needs to be transferred from the device to the cloud services.
- It can be done by passing MQTT messages to AWS IoT Services.
- It would require the setup of MQTT broker and client on RaspberryPi.

Wifi- SetUp

```
WiFiClientSecure espClient;
PubSubClient client(AWS_endpoint, 8883, callback, espClient); //set MQTT port number to 8883
#define BUFFER_LEN 256
long lastMsg = 0;
char msg[BUFFER_LEN];
int value = 0;
byte mac[6];
char mac_Id[18];

#define REPORTING_PERIOD_MS 1000

PulseOximeter pox;

void setup_wifi() {

    delay(10);
    // We start by connecting to a WiFi network

    espClient.setBufferSizes(512, 512);
    Serial.println();
    Serial.print("Connecting to ");
    Serial.println(ssid);

    WiFi.begin(ssid, password);

    while (WiFi.status() != WL_CONNECTED) {
        delay(500);
        Serial.print(".");
    }

    Serial.println("");
    Serial.println("WiFi connected");
    Serial.println("IP address: ");
    Serial.println(WiFi.localIP());

    timeClient.begin();
```

Reading sensor values

```
while 1==1:
    sleep(5)
    if connflag == True:
        ethName=getEthName()
        ethMAC=getMAC(ethName)
        macIdStr = ethMAC
        Ambient_Temp = sensor.get_ambient()
        Body_Temp = sensor.get_object_1()
        payloadmsg0="{\"
        payloadmsg1 = \"\\\"mac_Id\\\": \"\"
        payloadmsg2 = \"\\\", \\\"Ambient_Temp\\\": \"
        payloadmsg3 = \"\\\", \\\"Body_Temp\\\": \"\"
        payloadmsg4=}\"\"
        payloadmsg = \"{ } { } { } { } { } { } { }\".format(payloadmsg0, pa
        payloadmsg = json.dumps(payloadmsg)
        payloadmsg_json = json.loads(payloadmsg)
        mqttc.publish("IR_TEMPERATURE", payloadmsg_json , qos=1)
        print("msg sent: IR_TEMPERATURE" ) # Print sent temperature
        print(payloadmsg_json)
```

```
void loop() {

    if (!client.connected()) {
        reconnect();
    }
    client.loop();

    pox.update();

    long now = millis();
    if (now - lastMsg > 1000) {
        lastMsg = now;

        String macIdStr = mac_Id;
        float p = pox.getHeartRate(); // Get heart rate reading
        delay(10);
        float s = pox.getSpO2(); // Get spo2 reading
        snprintf (msg, BUFFER_LEN, \"{\\\"mac_Id\\\" : \\\"%s\\\", \\\"Pulse_R

        Serial.print("Publish message: ");
        Serial.println(msg);
        client.publish("Pulse_Oximeter", msg);
        Serial.print("Heap: "); Serial.println(ESP.getFreeHeap());
    }
}
```

MQTT

- Used for sending the data from devices to AWS IoT.
- Requires client and broker connection setup.
- The clients connected to the broker sends the messages.
- The broker distributes the messages to all the devices subscribed to it.

Setting Up the Client-Broker Connection

```
mqttc = paho.Client()
mqttc.connect()
mqttc.publish()
```

```
mqttc = paho.Client()
mqttc.on_connect = on_connect
mqttc.on_message = on_message
#mqttc.on_log = on_log

#### Change following parameters ####
awshost = "a2y9prqvlyko4a-ats.iot.us-east-2.amazonaws.com" # Endpoint
awsport = 8883 # Port no.
clientId = "Raspberrypi" # Thing_Name
thingName = "Raspberrypi" # Thing_Name
caPath = "/home/pi/AWSIoT/root-ca.pem"
certPath = "/home/pi/AWSIoT/certificate.pem.crt"
keyPath = "/home/pi/AWSIoT/private.pem.key" # <Th

mqttc.tls_set(caPath, certfile=certPath, keyfile=keyPath, cert_reqs=ssl.CER

mqttc.connect(awshost, awsport, keepalive=60) # connect to av

mqttc.loop_start() # Start the loo

while 1==1:
    sleep(5)
    if connflag == True:
        ethName=getEthName()
        ethMAC=getMAC(ethName)
        macIdStr = ethMAC
        Ambient_Temp = sensor.get_ambient()
        Body_Temp = sensor.get_object_1()
        payloadmsg0="{
        payloadmsg1 = "\"mac_Id\": \""
        payloadmsg2 = "\", \"Ambient_Temp\": \"
        payloadmsg3 = ", \"Body_Temp\": \""
        payloadmsg4="\"}"
        payloadmsg = "{} {} {} {} {} {} {} {}".format(payloadmsg0, payloadmsg1, payloadmsg2, payloadmsg3, payloadmsg4)
        payloadmsg = json.dumps(payloadmsg)
        payloadmsg_json = json.loads(payloadmsg)
        mqttc.publish("IR_TEMPERATURE", payloadmsg_json , qos=1) # to
        print("msg sent: IR_TEMPERATURE" ) # Print sent temperature msg on
        print(payloadmsg_json)

    else:
        print("waiting for connection...")
```

```

File ca = SPIFFS.open("/ca.der", "r"); //replace ca eith your uploaded file name
if (!ca) {
    Serial.println("Failed to open ca ");
}
else
    Serial.println("Success to open ca");

    delay(1000);

    if(espClient.loadCACert(ca))
        Serial.println("ca loaded");
    else
        Serial.println("ca failed");
        Serial.print("Heap: "); Serial.println(ESP.getFreeHeap());

    WiFi.macAddress(mac);
    snprintf(mac_Id, sizeof(mac_Id), "%02x:%02x:%02x:%02x:%02x:%02x",
    mac[0], mac[1], mac[2], mac[3], mac[4], mac[5]);
    Serial.print(mac_Id);
}

void loop() {

    if (!client.connected()) {
        reconnect();
    }
    client.loop();

    pox.update();

    Long now = millis();
    if (now - lastMsg > 1000) {
        lastMsg = now;

        String macIdStr = mac_Id;
        float p = pox.getHeartRate();        // Get heart rate reading
        delay(10);
        float s = pox.getSpO2();    // Get spo2 reading
        snprintf (msg, BUFFER_LEN, "{\"mac_Id\" : \"%s\", \"Pulse_Rate\" : %f, \"Spo2\" : %f}", macIdStr.c_str(), p, s);

        Serial.print("Publish message: ");
        Serial.println(msg);
        client.publish("Pulse_Oximeter", msg);
        Serial.print("Heap: "); Serial.println(ESP.getFreeHeap()); //Low heap can cause problems
    }
}

```

Challenges Faced

- Capturing the data simultaneously from two sensors.
- Setting up the connection between Rpi and AWS IoT.
- Displaying the data on a website.
- Setting up the wifi connection with the microcontrollers,
- Setting up the MQTT broker and Client connection.

Screenshots of the Webpage

Sensors Dashboard

Time-Stamp	Pulse_Rate	Spo2
11:57:57	59.642956	94
11:57:56	59.642956	94
11:57:55	59.642956	94
11:57:54	59.642956	94
11:57:53	59.642956	94
11:57:52	59.642956	94
11:57:51	59.642956	94
11:57:50	59.642956	94
11:57:49	59.642956	94
11:57:48	59.642956	94
11:57:47	59.642956	94
11:57:46	59.642956	94
11:57:45	59.642956	94
11:57:44	59.642956	94
11:57:43	59.642956	94
11:57:42	59.642956	94
11:57:41	59.642956	94
11:57:40	59.642956	94
11:57:39	59.642956	94
11:57:38	59.642956	94

